

# RC612 Multiphase Carbon and Hydrogen/Moisture Determinator

## Specification Sheet

### Instrument Range\* (250 mg for %; 7.2 in<sup>2</sup> strip sample for area)

% Carbon:	50 ppm to 20%
Area Carbon:	0.002 to 6.94 mg/in <sup>2</sup>
% H <sub>2</sub> O:	100 ppm to 20%

### Precision

Carbon:	25 ppm or 3% RSD of result; whichever is greater
H <sub>2</sub> O	50 ppm or 3% RSD of result; whichever is greater

### Calibration

Standards

### Detection Method

Solid-state infrared

### Sample Size

Powdered:	0.5 grams, nominal
Non-Powdered:	Strip or tubular samples up to 0.9 in. x 4 in. (23 mm x 102 mm)

### Chemical Reagents

- Sodium Hydroxide on Inert Base
- Anhydrous Magnesium Perchlorate
- Copper Oxide

### Gases (Carrier)

Oxygen, 99.5% pure, 40 psi (2.8 bar) ±10%  
Nitrogen, 99.5% pure, 40 psi (2.8 bar) ±10%

### Gas Flows

Oxygen:	3 lpm when in use
Nitrogen:	3 lpm when in use

### Furnace Range

Resistance Heated; near-ambient to 1100°C (2012°F)  
Programmable setpoints and ramp rates

### Environmental Conditions

Operating Temperature:	15°C to 35°C (59°F to 95°F)
Relative Humidity:	20% to 80%, non-condensing

### Electrical Power Requirements

230 V~ (±10%; at max load), 50/60 Hz, single phase, 30 A,  
23,600 BTU/hr

### Physical Dimensions\*\*

33 in. H x 22.5 in. W x 34.5 in. D (84 cm x 57 cm x 87 cm)

### Weight (Approximate)

255 lb. (116 kg)

### Part Numbers

RC612C	RC612 Determinator with Windows®-based software and external PC
RC612LC	RC612 Determinator with Windows-based software, external PC, and 50-sample Autoloader


### Options

620-400-210	Autoloader Add-On Kit
501-291	O <sub>2</sub> Regulator
764-216	Inert Regulator
751-300-100	L-250 Balance
619-995	Bar Code Reader (USB)
625-507	Kit Purifier

V~ denotes VAC.

\*Range may be affected by reducing or increasing sample mass; specified range based upon LECO standard materials and methods; range is matrix dependent.

\*\*Allow a 6-inch (15 cm) minimum access area around all units.



## Theory of Operation


The RC612 multiphase carbon and hydrogen/moisture determinator quantifies the carbon and hydrogen present in various organic and inorganic samples, and identifies the source of several types of carbon content.

The RC612 features a state-of-the-art furnace control system, which allows the temperature of the dual-stage furnace to be programmed from near-ambient to 1100°C.

Depending on the application, multiple furnace steps can be programmed by the operator and the furnace purged with oxygen or nitrogen, creating oxidizing or inert conditions in which the carbon and hydrogen present is combusted or reacted. An afterburner furnace (nominally set to 850°C) and a secondary oxidation catalyst are included in the flow path to ensure full combustion/reaction of all evolved components. Infrared detection is used to quantify the result either as a weight percentage or as a coating weight (mg/in<sup>2</sup>).

When combusted in an oxidizing atmosphere (O<sub>2</sub>), all forms of carbon (except some carbides such as SiC) are converted to CO<sub>2</sub>. In contrast, organic forms of carbon produce both H<sub>2</sub>O and CO<sub>2</sub>. Thus the presence of organic carbon may be verified by finding coincident peaks in H<sub>2</sub>O and CO<sub>2</sub>. Moisture and carbonate are detected when the sample is combusted in an inert (N<sub>2</sub>) atmosphere, with the furnace catalyst temperatures at 120°C. In this mode, organic carbon is normally not detected. Additional sources of carbon can often be differentiated by the temperature at which they oxidize or volatilize.

A slow ramping temperature program, from 100°C to 1000°C at 20°C per minute, can be used for the analysis of unknown samples. This type of analysis can be used to indicate the temperatures at which the differing forms of carbon are oxidized, thereby enabling the operator to optimize the furnace temperature program to provide more rapid quantitative results for each form of carbon present in this sample type.



Specifications and part numbers may change.  
Consult LECO for latest information.

3000 Lakeview Avenue • St. Joseph, MI 49085 • Phone: 800-292-6141 • Fax: 269-982-8977  
info@leco.com • www.leco.com • ISO-9001:2000 • No. FM 24045 • LECO is a registered trademark of LECO Corporation.

**World Headquarters/United States**  
In United States (800) 292-6141  
(269) 985-5496  
Outside U.S.A. (269) 983-5531

### International Offices/Subsidiaries

**Africa** • LECO Africa (Pty.) Ltd.  
Kempton Park, South Africa  
Ph: 27-11-974-1681 • Fax: 27-11-974-1848

**Argentina** • LECO Argentina S.A.  
Buenos Aires, Argentina  
Ph: 54-114-523-8077 • Fax: 54-114-523-8873

**Australia** • LECO Australia Pty. Ltd.  
Castle Hill, NSW, Australia  
Ph: 61-2-9894-5955 • Fax: 61-2-9894-5247

**Brazil** • LECO Instrumentos Ltda.  
Rio De Janeiro, Brazil  
Ph: 55-212-538-4250 • Fax: 55-212-538-4299

**Canada** • LECO Instruments Ltd.  
Mississauga, Ontario Canada  
Ph: 905-564-6577 • Fax: 905-564-6582

**China** • LECO Corp/Beijing Office  
Beijing, China  
Ph: 86-10-6512-8877 • Fax: 86-10-6528-0373

**Czech Republic** • LECO Instrumente Plzen s.r.o.  
Plzen/Czech Republic  
Ph: 420-37-751-0811 • Fax: 420-37-725-9304

**France** • LECO France  
GARGES-Les-GONESSE; Cedex, France  
Ph: 33-134-45-4600 • Fax: 33-139-86-4105

**Germany** • LECO Instrumente GmbH  
Mönchengladbach, Germany  
Ph: 49-2166-687-0 • Fax: 49-2166-687-100

**Hong Kong** • LECO Instruments Hong Kong Ltd.  
Kowloon, Hong Kong  
Ph: 852-2387-7028 • Fax: 852-2708-4388

**Italy** • LECO Italy, S.R.L.  
Milan, Italy  
Ph: 39-02-953-43391 • Fax: 39-02-953-43461

**Japan** • LECO Japan Corporation  
Tokyo, Japan  
Ph: 813 (5782) 7800 • Fax: 813 (5782) 7801

**Korea** • LECO Korea Co. Ltd.  
Gyeonggi-Do, Korea  
Ph: 82-31-478-2441 • Fax: 82-31-478-2440

**Malaysia** • LECO Instruments (M) Sdn. Bhd.  
Selangor, D.E. Malaysia  
Ph: 60-3-7805-2100 • Fax: 60-3-7805-2212

**Mexico** • LECO Mexico S.A. de C.V.  
Monterrey, NL CP 64900 Mexico  
Ph: 52-81-83-49-3927 • Fax: 52-81-83-49-6485

**The Netherlands** • LECO Europe B.V.  
Geleen, The Netherlands  
Ph: 31-46-4747473 • Fax: 31-46-4747333

**Poland** • LECO Polska Sp.zo.o  
Katowice, Poland  
Ph: 48-32-2000-760 • Fax: 48-32-2000-536

**Russia** • Z.A.O. LECO Center Moscow  
Moscow, Russia  
Ph: 7495-710-3818 • Fax: 7495-710-3826

**Spain** • LECO Instrumentos S.A.  
Madrid, Spain  
Ph: 3491-803-1250 • Fax: 3491-804-0577

**Sweden** • LECO Corporation Svenska AB  
Upplands Väsby, Sweden  
Ph: 468-594-11000 • Fax: 468-594-11011

**Taiwan** • LECO Instruments Taiwan Ltd.  
Taipei, Taiwan R.O.C.  
Ph: 886-22-518-4699 • Fax: 886-22-518-4671

**Thailand** • LECO Instruments (Thailand) Ltd.  
Bangkok, Thailand  
Ph: 66-2-886-4350 • Fax: 66-2-886-4351

**United Kingdom** • LECO Instruments (U.K.) Ltd.  
Stockport, Cheshire, England  
Ph: 44-161-487-5900 • Fax: 44-161-456-0969

**Venezuela** • Instrumentos LECO Venezuela, C.A.  
Puerto Ordaz, Edo. Bolivar, Venezuela  
Ph: 58-286-994-4835 • Fax: 58-286-994-1570

**Vietnam** • LECO (Vietnam) Co., Ltd.  
Ho Chi Minh City, Vietnam  
Ph: 848-829-8979 • Fax: 848-829-8978

**LECO Corporation**



Delivering the Right Results