

CHN628 Series Sulfur Add-On Module Specification Sheet



Instrument Range*	0.01 to 20 mg Sulfur	
Precision (Sulfur)	0.005 mg or 1% RSD (whichever is greater)	
Nominal Sample Weight	up to 350 mg, 250 mg nominal	
Analysis Time	60 to 120 seconds (typical)	
Detection Method	Infrared absorption	
Chemical Reagents	Magnesium Perchlorate (Anhydrous)	
Gas Requirements	Oxygen, 99.5% pure, 40 psi (2.8 bar) \pm 10% 3.5 lpm, 1.0 lpm purge	
Furnace	600 to 1450°C \pm 1% of setpoint; Horizontal resistance-type	
Operational Control	Windows®-based software for PC on CHN628	
Environmental Conditions	Operating Temp: 15°C to 30°C (59°F to 86°F) Humidity: 20% to 80%, non-condensing	
Physical Dimensions[‡]	31 in H x 13 in W x 27 in D (79 x 33 x 68 cm)	
Weight (approx.)	147 lb. (67 kg)	Shipping Weight (approx.): 198 lb. (543 kg)
Electrical Requirements	230 V~ (\pm 10%; at max load), 50/60 Hz, single phase, 25 A; 19,700 BTU/hr [†]	

Part Numbers

FP628SC	FP628 System for Nitrogen and Sulfur Determination; Includes PC, Monitor, Windows®-based software, module, and host system
FP628SLC	Same as FP628SC above, with Sulfur Module Autoloader
CN628SC	CN628 System for Carbon, Nitrogen, and Sulfur Determination; Includes PC, Monitor, Windows®-based software, module, and host system
CN628SLC	Same as CN628SC above, with Sulfur Module Autoloader
CHN628SC	CHN628 System for Carbon, Hydrogen, Nitrogen, and Sulfur Determination; Includes PC, Monitor, Windows®-based software, module, and host system
CHN628SLC	Same as CHN628SC above, with Sulfur Module Autoloader
628SADD	628 Series Sulfur Add-On Module for upgrading 628 Systems; Includes Windows®-based software, module, and necessary components for mating module to an existing system, (no PC or Monitor)
628SADDL	Same as 628SADD above, with Sulfur Module Autoloader

Optional Accessories

621-192	Dual Monitor Add-On Kit
528-203-250	Combustion Boats

*Adjusting sample size may extend instrument range.
[†]Average output based on nominal operating parameters.
[‡]Allow a 6-inch (15 cm) minimum access area around all units.
 V~ denotes VAC.


Theory of Operation

LECO offers a Sulfur Add-On Module for the CHN628 Series of Elemental Determinators. The add-on module provides sulfur analysis for any element combination of the CHN628 Series—FP (N), CN, or CHN. The 628 S module is specifically designed to determine the sulfur content in a wide variety of organic materials such as coal, coke, and fuel oils, as well as some inorganic materials such as soil, cement, and limestone.

Analysis begins as a sample is weighed into a combustion boat and placed in the furnace with pure oxygen typically regulated at 1350°C. Sulfur within the sample is evolved from the sample and forms SO₂. The sample gases exiting the furnace are first swept through the boat stop to the back of the inner combustion tube, then forward between the inner and outer combustion tubes. This allows the sample gases to remain in the high-temperature zone for a longer period and permit efficient oxidation. From the combustion system, the gases flow through an anhydrous tube to remove moisture, and through a flow controller that sets the flow of sample gases through the sulfur infrared detection cell within the CHN628 series instrument.

While the CHN628 Series instrument and sulfur add-on module can be loading/analyzing/operating completely independent of each other, the module requires the detection capabilities and PC offered within the CHN628 Series instrument system in order to complete the analysis.

Flow Diagram



Specifications and part numbers may change.
Consult LECO for latest information.

3000 Lakeview Avenue • St. Joseph, MI 49085 • Phone: 800-292-6141 • Fax: 269-982-8977
info@leco.com • www.leco.com • ISO-9001:2008 HQ-Q-994 • LECO is a registered trademark of LECO Corporation.

LECO Corporation

