

Supplies for Organic Analysis



LECO
EMPOWERING RESULTS

This page intentionally blank.







Table of Contents


Released May 4, 2021

828 Series.....	4
928 Series.....	7
TruSpec Micro.....	10
Micro Oxygen Add-On Module	13
SC832 Series	15
TGA801	17
TGM800	19
AC600	20
AC500	22
AF700.....	24
Ordering Information	26

REFERENCE MATERIALS AND STANDARDS

Current part numbers and values for LECO Reference Materials and Standards (209-296) can be found at <https://www.leco.com/consumables>







Consumables & Spare Parts

GENERAL COMBUSTION SUPPLIES

501-171-HAZ	Anhydron 454 g (1 lb.) bottle
502-973	Copper Stick, Deoxidized, 200 g
502-995	Copper Turnings, Degassed, 60 g
501-081	Glass Wool 454 g (1 lb.) box
763-265	Glass Wool 227 g (0.5 lb.) package
502-174-HAZ	Lecosorb 500 g bottle
502-049	N-Catalyst 50 g bottle
608-379	Quartz Wool Strips 10/pk
502-310	Steel Wool 454 g (1 lb.)
503-510-HAZ	Secondary Furnace Reagent, 100 g
622-002-319	Ballast C-Flex Replacement Assy (FP/CN models)
622-002-632	Upper Pinch Valve Tubing Assy (CHN model)
622-002-633	Lower Pinch Valve Tubing Assy (CHN model)

CAPSULES AND CRUCIBLES

502-338	Gel Capsule, Small 400/pkg
502-382	Gel Capsule, Medium 400/pkg
502-810	Gel Capsule, Large 400/pkg
614-961-110	Porous Crucibles for solids 10/pk
617-605	Porous Crucibles for liquids 10/pk
502-040	Tin Capsules, Small 100/bottle
502-040-100	Tin Capsules, Small 1000/bottle
502-167	Tin Capsules, Medium 100/bottle
502-825	Tin Capsules, Large 50/bottle
502-186-200	Tin Foil Cups, Small 100/bottle
502-186-100	Tin Foil Cups, Small 1000/bottle
502-397	Tin Foil Cups, Large 100/bottle
502-397-400	Tin Foil Cups, Large 400/bottle

OPERATION ACCESSORIES

502-007	Tube Cleaner
601-981	Wire Brush
616-152	Crucible Extractor Tool
775-306	Secondary Filter, 10 micron
619-591-699	Filter In-Line, 0.1 micron
502-023	Funnel Quick Disc
778-321	Lance Extraction Tool
502-213	Leak Detection Solution
762-515	O-ring Tool - Removal
616-513	O-ring Tool - Installation
604-398	Sample Cup Holder, Large
604-373	Sample Cup Holder, Small
778-405	10-Sample Storage Block
501-614	Spatula
619-180-110	Stackable Carousel Assy
625-603-404	Tube Packing Tool
502-028	Tweezers
604-378	Tweezers, Curved (SST)
501-241	Vacuum Grease

COMBUSTION TUBES AND GLASSWARE

619-316	Glass Filter Tube (wide mouth top, small bottom) requires one (1) 765-976 Top O-ring and one (1) 772-520 Bottom O-ring for sealing; (FP/CN models)
619-425	Straight Filter Glass Tube (CHN model)
780-899	Honeycomb Ceramic Stop
619-065	Combustion Tube
619-591-740	Aliquot Scrubber Tube
625-603-030	Reduction Tube
616-146	Lance Tube
625-603-353	Incoming TC Carrier Gas Scrubber Tube
622-002-748	Top Load Secondary Furnace Reagent Tube (CHN model)
619-777	Middle Directing Tube (CHN model)

MAINTENANCE PARTS KITS

104-196	CN/FP828 Maintenance Parts Kit
104-206	CHN828 Maintenance Parts Kit

5K CONSUMABLE PACKS


502-956-HAZ	828 Series Small Tin Foil Cup
502-957-HAZ	828 Series Large Tin Foil Cup
502-958-HAZ	828 Series Medium Gel Cap
502-959-HAZ	828 Series Small Gel Cap
502-960-HAZ	828 Series Small Tin Cap











1K CONSUMABLE PACKS

502-951-HAZ	828 Series Small Tin Foil Cup
502-952-HAZ	828 Series Large Tin Foil Cup
502-953-HAZ	828 Series Small Gel Cap
502-954-HAZ	828 Series Medium Gel Cap
502-955-HAZ	828 Series Small Tin Cap

REFERENCE MATERIALS AND STANDARDS

Current part numbers and values for LECO Reference Materials and Standards (209-296) can be found at <https://www.leco.com/consumables>







Consumables & Spare Parts

GENERAL COMBUSTION SUPPLIES

501-171-HAZ	Anhydron 454 g (1 lb) bottle
502-973	Copper Stick Deoxidized, 200 g
502-995	Copper turnings Degassed, 60 g
501-081	Glass Wool 454 g (1 lb) box
626-001-871	Furnace Reagent Tube (primary/secondary filter tube) (O-ring top 773-320, O-ring bottom 772-520)
625-603-030	Reduction Tube (O-ring top 773-320, O-ring bottom 772-520)
775-306	Screen Filter
619-591-699	Filter In-Line 0.01 Micron Disposable
502-174-HAZ	Lecosorb 500 g bottle
502-049	N-Catalyst 50 g bottle
502-310	Steel Wool 454 g (1 lb)
606-308	Boat Stop, Ceramic
501-426	ComAid™ for Solids 454 g (1 lb) bottle
501-427	ComAid™ for Liquids 454 g (1 lb) bottle
528-203-250	Combustion Boat - Ceramic 250/pkg
625-505-430	Nickel Boat (50/pkg)
502-343	Nickel Boat Liner 100/pkg
625-603-331	Pinch Valve Replacement Tubing Assembly
601-691	C-Flex Tubing (25 ft)
625-602-297	O-ring for Aliquot Doser Stem
611-476	O-ring Inner Combustion Tube
611-477	O-ring Outer Combustion Tube
772-520	O-ring Lance
625-601-515	Combustion Tube Insulation
625-601-554	Outer Combustion Tube
625-601-555	Inner Combustion Tube
625-602-187	Lance Tube
625-525-129	Heating Element Assembly w/Clamp and Lugs

OPERATION ACCESSORIES

761-929	Crucible Tongs
626-002-644	Heating Element Connector Tool
502-007	Tube Cleaner
625-603-399	Rectangular Foam Swab 100/pk
501-241	Vacuum Grease
502-213	Leak Detect Solution
501-614	Spatula
501-062	Boat Puller
502-023	Funnel
616-513	O-ring Installation Tool for Aliquot Doser Stem
606-403	Tube Brush

MAINTENANCE PARTS KITS

104-194	CN/FP928 Maintenance Parts Kit
104-203	CNS928 Maintenance Parts Kit

5K CONSUMABLE PACKS

502-726-HAZ	FP/CN928 Solid
502-727-HAZ	FP/CN928 Liquid
502-729-HAZ	CNS/NS928


1K CONSUMABLE PACKS

502-724-HAZ	FP/CN928 Solid
502-725-HAZ	FP/CN928 Liquid
502-728-HAZ	CNS/NS928





REFERENCE MATERIALS AND STANDARDS

Current part numbers and values for LECO Reference Materials and Standards (209-296) can be found at




<https://www.leco.com/consumables>





TruSpec Micro



ALTERNATIVE PRIMARY FURNACE PACKING







Consumables & Spare Parts

GENERAL COMBUSTION SUPPLIES

501-171-HAZ	Anhydron 454 g (1 lb.) bottle
502-190	Copper Oxide Wire 100 g bottle
502-705	Copper Sticks 100 g ampule
502-878	Copper Sticks 100 g bottle
502-995	Copper Turnings 60 g bottle
501-081	Glass Wool 454 g (1 lb.) box
763-265	Glass Wool 227 g (0.5 lb.) package
502-174-HAZ	Lecosorb 500 g bottle
502-644	Micro Furnace Reagent
502-049	N-Catalyst 50 g bottle
608-379	Quartz Wool Strips 10/pk
502-301	Tungsten Trioxide (125 g)

CAPSULES AND CRUCIBLES

614-961-110	Porous Crucibles for solids 10/pk (use only with 614-962 lance)
502-293	Silver Capsules for Micro Analysis 500/bottle
502-206	Silver Capsules for Micro Analysis 1000/bottle
502-227	Tin Capsules for Micro Analysis 500/bottle
502-489	Reticulated Micro Crucibles 10/pk (for CHNS only)
625-511-154	Porous Reticulated Crucibles (small) 5/pk

OPERATION ACCESSORIES









502-005	Brush
601-981	Wire Brush
616-152	Crucible Extractor Tool
775-306	Secondary Filter, 10 micron
775-307	Filter Tool, Thumbscrew
502-023	Funnel Quick Disc
778-321	Lance Extraction Tool
502-213	Leak Detection Solution
603-350	Tool O-ring Installation
762-515	Tool O-ring Removal
604-398	Sample Cup Holder, Large
604-373	Sample Cup Holder, Small
778-405	10-Sample Storage Block
501-614	Spatula
502-007	Tube Cleaner
502-028	Tweezers
604-378	Tweezers, Curved
501-241	Vacuum Grease
619-380	Vacuum Cleaner Attachment Assembly

COMBUSTION TUBES AND GLASSWARE

619-425	Micro In-Line Anhydron Tube
774-638	Reagent Tube, requires 772-520 O-rings for sealing
775-601	Reagent Tube, requires 772-520 O-rings for sealing
619-154	Reduction Tube, requires 601-442 O-rings for sealing

REFERENCE MATERIALS AND STANDARDS

Current part numbers and values for LECO Reference Materials and Standards (209-296) can be found at <https://www.leco.com/consumables>





Consumables & Spare Parts

GENERAL COMBUSTION SUPPLIES

501-171-HAZ	Anhydron 454 g (1 lb.) bottle
501-170	Rare Earth Copper Oxide
502-196	Carbon Black 10 g bottle
502-878	Copper Sticks 100 g bottle
502-995	Copper Turnings 60 g bottle
501-081	Glass Wool 454 g (1 lb.) box
763-265	Glass Wool 227 g (0.5 lb.) package
606-026	Graphite Crucible 10/package
502-228	Graphite Pad 1 each
502-174-HAZ	Lecosorb 500 g bottle
601-963	Tin Capsule 1000/bottle (3.5 x 4 mm)
502-227	Tin Capsule 500/bottle (3.2 x 4 mm)

OPERATION ACCESSORIES

502-005	Brush
601-981	Wire Brush
775-306	Secondary Filter, 10 micron
502-213	Leak Detection Solution
603-350	Tool O-ring Installation
762-515	Tool O-ring Removal
608-135	Pyrolysis Tube Extractor Tool Assembly
604-398	Sample Cup Holder, Large
604-373	Sample Cup Holder, Small
778-405	10-Sample Storage Block
502-007	Tube Cleaner
502-028	Tweezers
604-378	Tweezers, Curved
501-241	Vacuum Grease
601-964	Micro Pliers
619-726	Micro Forceps
619-725	Micro Spatula

COMBUSTION TUBES AND GLASSWARE

608-065	Graphite Lance
608-130	Inner Pyrolysis Tube, requires one (1) 608-089 O-ring
619-154	Oxygen Catalyst Heater Tube
625-401-360	Outer Pyrolysis Tube, requires one 608-090 O-ring, and one 608-091 O-ring
775-601	Reagent Tube, requires 772-520 O-rings for sealing

REFERENCE MATERIALS AND STANDARDS

Current part numbers and values for LECO Reference Materials and Standards (209-296) can be found at <https://www.leco.com/consumables>



GENERAL COMBUSTION SUPPLIES

501-171-HAZ	Anhydrene 10-16 Mesh 1 lb.
501-426	COM-AID Solids and Powders 1 lb.
501-427	Com-Aid Liquids 1 lb.
502-321	Accelerator COM-CAT 250 g
529-203	Boat Combustion Zp P 100/package
528-203-250	Boat Crucible Zir 250/package
501-081	Wool Glass 1 lb.
763-265	Wool Glass 0.50 lb.
502-301	Tungsten Trioxide 125 g
769-608	Metal Trap Antimony 200 g, Used in Halogen Scrubber Packing
769-610	F-CL Absrb -8 +20 mesh 225 g, Used in Halogen Scrubber Packing

**COMBUSTION TUBES, GLASSWARE,
and HEATING ELEMENT**

606-308	Stop Ceramic Stepped
626-001-691	Tube Lance Manual Load, Requires 625-401-302 O-Ring
626-001-692	Tube Purge Sst Manual, Requires 625-401-301 O-Ring
625-601-554	Tube Combustion Outer, Requires 611-477 O-Ring
625-601-555	Tube Combustion Inner, Requires 611-476 O-Ring
626-002-572	High Temperature Outer Combustion Tube
626-002-573	High Temperature Inner Combustion Tube
625-525-129	Heating Element Assy W/Clamp and Lugs
626-001-871	Tube Boros Reagent 1.38 OD 7.0
626-001-750	Pack O-Ring Spare SC/C/S832 Contains: (8) 772-520 O-Ring 108 (.250x .436x.093v); (4) 773-320 O-Ring 121 1.062x 1.248x.093v; (2) 611-476 O-Ring 326 (1.625x 1.999x.187a); (2) 611-477 O-Ring 330 (2.125x 2.499x.187a); (2) 625-401-301 O-Ring 006 (.114x .184x .070a); (2) 625-401-302 O-Ring 107 (.206x .309x.103a)

OPERATION ACCESSORIES

501-062	Assy Boat Puller
761-929	Tong Crucible
501-614	Spatula Sample
503-032	Scoop Glass Accelerator Small
619-591-699	Filter In-Line 0.1 Micron Disposable*
502-007	Cleaner Tube 8 mm x 4 ft
501-241	Lub Grease Vacuum 5.3 oz
502-023	Funnel Quick Disc Tube
776-978	Scoop Glass 1.5 ml, Large
633-103-356	Kit Regulator Oxygen CGA 540
751-350-110	Kit Balance 4-plc Secura124-1s
611-844	Spatula Flat Spoon Sst
626-002-251	Kit Pump Rebuild

MAINTENANCE PARTS KITS

104-186	832 Series w/Autoloader
104-188	832HT Series
104-190	832HT Series w/Autoloader


5K CONSUMABLE PACKS

502-771-HAZ	832 Series
-------------	------------

REFERENCE MATERIALS AND STANDARDS

Current part numbers and values for LECO Reference Materials and Standards (209-296) can be found at
<https://www.leco.com/consumables>


LECO Consumables & Spare Parts



TGA801



Consumables & Spare Parts



OPTIONAL ACCESSORIES


619-593-570 Scanner Bar Code Laser Bluetooth
619-995 Scanner Bar Code Laser USB
621-011-553 Weight Calibration 50 g

CRUCIBLES AND SPARE PARTS

621-331 Ceramic Crucible (5 pkg)
529-048 Ceramic Crucible Cover (5 pkg)
502-628 Crucible Foil Liner (1200/pkg)
621-010-236 Crucible Foil Alum 1.51 Dia 1000/pkg
621-011-378 Crucible Holder Tray
778-898 Crucible Cover Forceps
621-223 Filter Air .38t x 3.7 Dia.
621-011-549 Filter Air Foam 3.88 x 7.86 x .5 45 ppi Gry
619-591-699 Filter In-Line 0.1 Micron Disposable
621-011-507 Double Ended Scoop
614-442-110 Turntable Crucible Large Annealed

REFERENCE MATERIALS AND STANDARDS

Current part numbers and values for LECO Reference Materials and Standards (209-296) can be found at <https://www.leco.com/consumables>



OPTIONAL ACCESSORIES


621-011-553	Weight Calibration 50 g
619-995	Scanner Bar Code Laser USB
619-593-570	Scanner Bar Code Laser Bluetooth
621-011-560	Kit Oven Temp Verification
621-011-574	Kit Extension Exhaust TGM

CRUCIBLES AND SPARE PARTS

502-237	Wipes Foam 25/pk
621-011-240	Filter Air Foam 6.00 x 11.00 x 1.00 25 ppi
621-011-241	Filter Air Foam 6.00 x 11.00 x 1.00 10 ppi
621-011-103	Filter Air Foam 6.00 x 11.00 x 1.00 45 ppi
621-010-609-110	Assy Balance Complete TGM Tested
621-010-236	Crucible Foil Alum 1.51 dia 1000/pk
621-010-956	Crucible Foil Alum 2.40 dia 1000/pk

REFERENCE MATERIALS AND STANDARDS

Current part numbers and values for LECO Reference Materials and Standards (209-296) can be found at <https://www.leco.com/consumables>



GENERAL COMBUSTION SUPPLIES

774-649	Cleaning Brush Assembly
774-204	Crucible 10/pkg.
502-461	Cotton Thread Fuse, 10 cm, 375/pkg
502-462	Fuse Wire, 10 cm, 375/pkg
621-605-803	Fixed Wire Igniter Assembly, 5/pkg
502-460	Fixed Wire Igniter

OPERATION ACCESSORIES

501-241	Vacuum Grease (Silicone)
633-101-923	AC600 Vessel Holder
621-657	Filter Cartridge, 1 micron
621-605-744	Flush Adapter
762-515	O-Ring Removal Tool
776-978	Glass Scoop
502-051	Tweezers, Curved Fine Point
621-605-517	Vessel Lid Removal Tool
621-605-750	Vessel Tool Kit

OPTIONS





751-350-110	Electronic Balance
621-322	Pellet Press
621-606-160	External Chiller Kit, 120V, 60 Hz (sized for single AC600)
621-606-162	External Chiller Kit, 240V, 50 Hz (sized for single AC600)
621-606-164	External Chiller Kit, 120V, 60 Hz (sized for dual AC600)
621-606-166	External Chiller Kit, 240V, 50 Hz (sized for dual AC600)
633-103-354	Compressed Air Regulator
619-995	Bar Code Scanner USB
621-453-110	Deskjet Printer Kit
621-605-958	AC600 Prep Station
621-605-686	O-Ring Pack for Vessel
621-605-687	O-Ring Pack for Prep Station
621-605-979	Halogen Resistant Combustion Vessel II
621-605-978	Standard Combustion Vessel II

O-RINGS

611-543	Prep Station Nipple Charge
621-605-560	Vessel Cap
621-605-664	Check Valve
621-605-973	Seal Insert Exhaust
772-738	Prep Station Valve Inserts
773-913	Electrode Post
776-562	Ceramic Washer
621-656-001	Water Filter Housing

REFERENCE MATERIALS AND STANDARDS

Current part numbers and values for LECO Reference Materials and Standards (209-296) can be found at <https://www.leco.com/consumables>





Consumables & Spare Parts

GENERAL COMBUSTION SUPPLIES


502-892	Benzoic Acid Tablets 50/pkg.
774-649	Cleaning Brush Assembly
502-682	Coal Prox Plus High Standard Ref. Matl.
774-204	Crucible 10/pkg.
502-296	Filters 100/pkg.
776-465	Air Filter
608-671	Particle Filter
776-978	Glass Scoop
502-266	Fuse Wire
502-461	Cotton Thread Fuse, 10 cm, 375/pkg.
502-460	Fixed Wire Igniter, 5/pkg.
502-462	Fuse Wire, 10 cm, 375/pkg.
621-245	Vessel Thread Igniter Kit
621-246	Vessel Wire Igniter Kit
608-778	Fuse (3.15A, 115V)
608-779	Fuse (1.6A, 230V)
608-734	Combustion Vessel Fuse Lead (Combustion Vessel Harness Assembly)
608-761-002	Bottom Outlet Pipette Tubing
608-762-004	Bottom Pipette Tubing
608-657	5 Feet Replacement Tube
608-749	Screen (Filter)
608-751	Filter Retaining Ring
501-241	Vacuum Grease (Silicone)
502-121	Wetting Agent

OPTIONS

751-350-110	Electronic Balance
608-948	Combustion Vessel Cleaning Station
621-322	Pellet Press
601-483-100	Paper (2400 sheets)
601-480	Dot Matrix Printer

REFERENCE MATERIALS AND STANDARDS

Current part numbers and values for LECO Reference Materials and Standards (209-296) can be found at <https://www.leco.com/consumables>





Consumables & Spare Parts

GENERAL OPERATING SUPPLIES

625-501-431	Door Interlock Bypass Key
761-929	Tongs
625-501-662	Ash Cone Tray (Dual sided, 5 and 6 place, 100/pkg)
776-957	Ash Cone Mold, ASTM D1857
625-501-363	Ash Cone Tray Tool
502-496	Gold Wire Calibration Standard (6 in. x 0.2 in. diameter)
501-951	Nickel Wire Calibration Standard (6 in. x 0.2 in. diameter)
502-010	Dextrin Solution
501-233-HAZ	Ceramic Alumina Adhesive
625-501-471	Test Pin Tray (Metal)
763-265	Glass Wool (0.5 lb)
616-504	Filter (0.01 micron)

FURNACE AREA AND GLASSWARE

613-564	T-Clamp Expansion Tool
625-501-328	AF700 Replacement Heating Element Kit including: (2) 633-101-361 Braided Wire w/Clip (5) 633-101-362 Braided Wire (2) 625-501-149 O-Ring (12) 633-103-092 Ceramic Fiber Washer (3) 625-501-430 O-Ring (1) 625-501-464 Heating Element (3/4 inch—6 in set) (12) 625-501-470 Heating Element Clamp (3/4 inch)
625-501-361	AF700 Replacement Furnace Tube Kit including: (2) 625-501-149 O-Ring (1) 625-501-161 O-Ring (2) 625-501-382 Ring Insulation (1) 625-501-391 Mullite Furnace Tube (3) 625-501-430 O-Ring
633-103-027	AF700 High Temperature Kit including: (12) 190-728 Screws (12) 193-055 Washer (2) 625-501-149 O-Ring (12) 633-103-092 Ceramic Fiber Washer (2) 625-501-382 Ring Insulation (3) 625-501-430 O-Ring (1) 625-501-464 Heating Element Set (3/4-inch—6 in set) (2) 633-101-991 SST Braided Wire w/Clip (5) 633-101-992 SST Braided Wire (12) 194-172 Hex Nut (1) 633-101-997 Alumina Furnace Tube
625-501-394	D-Tube
633-101-990	Mullite Furnace Gas Inlet
633-101-441	Furnace Thermocouple

REFERENCE MATERIALS AND STANDARDS

Current part numbers and values for LECO Reference Materials and Standards (209-296) can be found at <https://www.leco.com/consumables>

Ordering Information/General Terms

Ordering, Customer Assistance, and all Other Inquiries



Order Online NOTE: Available to USA customers only.

- Visit www.lecodirect.com to register for an account



Phone

800-292-6141 or
269-985-5496



Email

sales@leco.com
for purchase orders or
pricing information

info@leco.com
for general inquiries



Fax

269-982-8987



Mail

LECO Corporation
3000 Lakeview Avenue
St. Joseph, MI 49085

VISA, Mastercard, and American Express Accepted.



Terms and Conditions

- Terms: Net 30 days with approved credit.
- Prices shown are US dollars, F.O.B. St. Joseph, Michigan. Freight charges are prepaid and added to the invoice.
- Sales taxes will be added where applicable.
- Returns must be approved in advance; re-stocking fee may be applied. Return Material Authorization (RMA) is required for all returns. Contact LECO for RMA information.
- Disclaimer: LECO products are intended for use by professional laboratory personnel who are familiar with the handling of these materials. Our products are not for use in household, drug, medical, or cosmetic applications. The buyer and/or user assumes full responsibility for the safe handling, storage, application, and proper disposal of all products ordered from this catalog.
- Complete terms and conditions are included with quotations and invoices.
- **Prices, part numbers, and availability are subject to change without notice.**

This page intentionally blank.



3000 Lakeview Avenue | St. Joseph, MI 49085 | 800-292-6141 | Phone: 269-985-5496
info@leco.com | www.leco.com | ISO-9001:2015 Q-994 | LECO is a registered trademark of LECO Corporation.

LECO Corporation

LECO

EMPOWERING RESULTS

© 2021 LECO Corporation